# CytoBrite® Slide Incubation System USER MANUAL

Cat. #2019-00-1 (115/230V)

**FOR RESEARCH USE ONLY**

---

## Warranty

SciGene warrants that the CytoBrite® Slide Incubation System described in this manual shall be free of defects in materials and workmanship for a period of 12 months from date of delivery. This warranty does not cover removable slide trays or accessories. In the event of a defect during the warranty period, SciGene's limit of liability will be to, at its sole discretion: (a) provide replacement parts; (b) perform factory repairs; or (c) replace the product. The foregoing warranty is void in the event the unit was abused or modified or used in a manner inconsistent with its intended purpose. SciGene makes no other warranty, expressed or implied including warranties of merchantability and fitness for a particular purpose. In no event shall SciGene be liable for any direct, indirect, special, incidental or consequential damages or for any damages resulting from loss arising out of or in connection with the sale, use or performance of the product.

## Copyright

Copyright ©2014-2016 SciGene Corporation. All rights reserved. SciGene, CytoBrite and CytoBond are trademarks of SciGene Corporation, Sunnyvale, CA. All other trademarks used in this manual are the property of their respective owners.

## Table of Contents

|  |           |
|--|-----------|
| <b>I. SAFETY NOTICES</b> .....                             | <b>3</b>  |
| A. Intended Use.....                                       | 3         |
| B. Instrument Safety .....                                 | 3         |
| C. Symbols and Conventions.....                            | 3         |
| D. Warnings .....  | 4         |
| E. Cautions .....  | 5         |
| F. Compliance .....  | 5         |
| <b>II. UNPACKING AND SET UP</b> .....                      | <b>7</b>  |
| A. Unpacking the Instrument .....                          | 7         |
| B. Parts Provided.....                                     | 7         |
| C. Coverslip Sealant Provided .....                        | 7         |
| D. Environmental Requirements.....                         | 7         |
| E. Installing the Instrument .....                         | 7         |
| <b>III. USING YOUR CYTOBRITE SYSTEM</b> .....              | <b>8</b>  |
| A. System Components .....                                 | 8         |
| B. Powering On / System Preheat.....                       | 9         |
| C. Preparing and Loading Slides / Starting a Protocol..... | 9         |
| D. Using the Controller .....                              | 9         |
| E. Using the USB Drive .....                               | 9         |
| <b>IV. CREATING AND EDITING PROTOCOLS</b> .....            | <b>10</b> |
| A. Creating a Protocol .....                               | 10        |
| B. Editing a Protocol.....                                 | 11        |
| <b>V. ONETEMP MODE</b> .....                               | <b>12</b> |
| <b>VI. RUNNING A PROTOCOL</b> .....                        | <b>13</b> |
| <b>VII. MANAGING FILES</b> .....                           | <b>14</b> |
| A. Viewing a Protocol .....                                | 14        |
| B. Copying a Protocol.....                                 | 14        |
| C. Renaming a Protocol.....                                | 15        |
| D. Deleting a Protocol .....                               | 15        |
| <b>VIII. SETTINGS AND PREFERENCES</b> .....                | <b>16</b> |
| A. Datalogging to USB Key.....                             | 16        |
| B. Calibrating Slide Block Temperature .....               | 17        |
| C. Upgrading Software .....                                | 18        |
| <b>IX. MAINTENANCE</b> .....                               | <b>19</b> |
| A. Cleaning the Instrument .....                           | 19        |
| B. Cleaning Air Vents.....                                 | 19        |
| C. Checking and Replacing Fuses .....                      | 19        |
| <b>X. TROUBLESHOOTING</b> .....                            | <b>20</b> |
| <b>XI. SPECIFICATIONS</b> .....                            | <b>20</b> |
| <b>XII. ORDERING INFORMATION</b> .....                     | <b>21</b> |
| <b>XIII. DECLARATION OF CONFORMITY</b> .....               | <b>22</b> |

## I. SAFETY NOTICES

### A. Intended Use

The **CytoBrite® Slide Incubation System** is intended for the heating and cooling of biological samples. The instrument should only be used according to instructions provided in this User Manual and other SciGene technical documents. If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.

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### B. Instrument Safety

Before operating the instrument, read the information in this section concerning hazards and potential hazards. Ensure that anyone involved with the instrument's operation is instructed in both general safety practices for laboratories and specific safety practices for the instrument.

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### C. Symbols and Conventions

The following chart is an illustrated glossary of the electrical symbols used on the **CytoBrite System**. Whenever such symbols appear on instruments, please observe appropriate safety measures.

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#### 1. Electrical Symbols

This symbol indicates that this is a protected ground terminal that must be connected to earth ground before any other electrical connections are made to the instrument.

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CAUTION: This symbol alerts you to consult this Operator's Manual for further information and to proceed with caution.

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This symbol indicates the OFF position of the main POWER switch.

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This symbol indicates the ON position of the main POWER switch.

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## 2. Non-Electrical Symbols

CAUTION: This symbol illustrates a heat hazard. Proceed with caution when working around these areas to avoid being burned by hot components.

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## D. Warnings

Failure to comply with the following warnings that are affixed to the product can lead to possible personal injury or death.

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This symbol on the rear of the instrument indicates the presence of the fuse box. **Warning: For Continued Protection Against Fire, Replace Only with Same Type Rating of Fuse.** Always disconnect the power cord before attempting to replace the fuse.

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## E. Cautions

Failure to comply with the following cautionary statement affixed to the product may lead to possible personal injury.

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This symbol located both on and in front of the lid indicates the potential presence of a Hot Surface. Use care when working in this area to avoid being burned.

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This symbol located on the back of the instrument warns the user to select the correct voltage before connecting the power cord. Operating with the voltage selector switch in the wrong position will damage the instrument and create a safety hazard.

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## F. Compliance

### 1. European Conformity (CE)

This symbol indicates the instrument is in compliance with all applicable European Union Electromagnetic Compatibility and Low Voltage Directives.

### 2. Waste Electrical & Electronic Equipment Directive (WEEE)

This symbol indicates the instrument is required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive 2002/96/EC, providing environmentally safe disposal of end of life equipment through recycling.

Contact your institution to provide environmentally safe disposal. If this is not available, contact your local SciGene distributor or SciGene Technical Support [techserv@scigene.com] for assistance. Do not treat electrical and electronic equipment as unsorted municipal waste.

### RoHS 3. Restriction of Hazardous Substances (RoHS)

This symbol indicates the instrument is in compliance with the European Union's (RoHS 2) Directive 2011/65/EU on the restriction of the use of certain hazardous substances in electrical and electronic equipment (EEE).

#### 4. Electromagnetic Compatibility (EMC)

The CytoBrite System is a Class A digital device under FCC Title 47 Part 15B and designated as Class A electrical equipment for measurement, control, and laboratory use (EN61326).

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#### 5. FCC Warnings and Notes

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. The instrument qualifies as an “exempted device” under 47 CFR 15.103(c), in regard to the cited FCC regulations in effect at the time of manufacture. Shielded cables must be used with this unit to ensure compliance with the Class A FCC limits.

## II. UNPACKING AND SET UP

### A. Unpacking the Instrument

To unpack the instrument, open the shipping carton and remove the two foam inserts. Take out the bag containing the manual, power cord and USB key. Remove the instrument from the box by reaching down both sides to lift from the base and place it on the bench. Open the lid to verify two slide trays are present. Inspect the instrument and accessories. *If any damage is evident, retain all shipping materials and contact SciGene [custserv@scigene.com] or your local distributor.*

### B. Items Provided

- Heating/cooling unit
- 2x Slide Trays
- Power Cord
- User Manual
- USB Key

### C. Coverslip Sealant Provided

One complimentary 100 ml bottle of CytoBond® Removable Coverslip Sealant (SciGene cat. # 2020-00-1) is shipped separately from the instrument. Use of CytoBond Sealant is required for sealing coverslips; to prevent probe evaporation without humidification. Contact SciGene [custserv@scigene.com] or your local distributor to order additional bottles.


### D. Environmental Requirements

Ensure that the area where the **CytoBrite System** is installed meets the following conditions, for reasons of safety and performance:

|                     |  |
|---------------------|--|
| Ambient temperature | +15 to +32°C (58 to 90°F)              |
| Relative humidity   | 20 to 80% RH non-condensing            |
| Air flow clearance  | 6 inches (15 cm) minimum on both sides |

### E. Installing the Instrument

Place the instrument on a level surface within a few feet of the power source. Position the system so the lid can be easily opened without interference. Allow 18 inches of vertical clearance and 16 inches of horizontal clearance. Ensure 3 inches of clearance at the back to attach the power cord. Using only the power cord provided, plug it into the back of the unit and then to a properly grounded outlet. Turn on power using the ON/OFF switch on the front control panel.


## III. USING YOUR CYTOBRITE SYSTEM


### A. System Components

| Component                            | Function  |
|--------------------------------------|---|
| SlideSense™ Temperature Sensor       | Senses and transmits slide block temperature information to controller              |
| Temperature Controller               | Used to program and monitor instrument operation and report slide block temperature |
| Slide Tray                           | Holds slides during sample preparation and incubation                               |
| Power Switch                         | Turns on main power to unit   |
| USB Port                             | Used for saving datalog files and installing software                               |
| Hinged Lid                           | Covers slide trays during processing  |
| Block Temperature Calibration Jack   | Connection for digital thermometer cable to calibrate slide block temperature       |
| Power Input Module and Fuse Housing  | Location of main fuses and power cord connection                                    |
| CytoBond Removable Coverslip Sealant | Seals coverslips to slides  |

### C. Powering On / System Preheat

To turn ON the instrument, press the top of the power switch to the “I” position. The display briefly shows the software version number then automatically starts a five minute 37C pre-heating step. When the preheat step is completed, the main menu appears with the slide block temperature in the upper right corner. The instrument is now ready for use.

### D. Preparing and Loading Slides / Starting a Protocol

Slides are prepared for molecular analysis following your standard procedures using the CytoBrite Slide Trays provided. CytoBond® Removable Coverslip Sealant (SciGene cat. # 2020-00-1) is used in place of rubber cement for sealing coverslips to prevent probe evaporation without humidification.

1. Place slides to be analyzed in the tray(s), pipet probes and place coverslips.
2. Seal coverslips by dispensing a continuous bead of CytoBond Sealant around edges.
3. Open the instrument lid and insert the trays(s)
4. Select and start the desired program (See section V for how to create programs.)
5. The instrument heats the slides to the programmed temperatures and times then returns to 37C until the program is ended by the operator.


### E. Using the Controller

A 4-line LED display controller is used to CREATE, RUN, EDIT and VIEW protocols, manage FILES (rename, copy, delete) and customize SETTINGS (including calibrating temperature and choosing alert sounds).

Navigating menus and selecting options is performed with four buttons on the display panel:

- **Home** — returns to the main menu screen.
- **Select** — chooses the option highlighted in the display.
- **Forward** — increments or toggles between options.
- **Back** — returns to a previous screen or option.

Details on using controller menus can be found in later sections of this manual.

### F. Using the USB Drive

A USB key inserted into the USB DRIVE can be used to:

- Capture a time and temperature log when running a protocol
- Allow protocols to run directly from a USB key.
- Provide additional capacity for storing protocols beyond the 10 file storage limit on the controller.
- Upload new controller software.


## IV. CREATING AND EDITING PROTOCOLS

### A. Creating a Protocol

Protocols are created and saved in memory on the controller (10 files maximum) or a USB key. Names are 8 alphanumeric characters in length including spaces.

1. From the Home screen, select Create.
2. Enter the name of the protocol by scrolling through the alphanumeric characters (space: A-Z; 0-9) using the Forward button. Sequentially select a desired character using the Select button. Names must be 8 characters in length, including spaces.
3. On the following screens enter the desired denature temperature then the time in minutes and seconds.
4. Save the completed protocol to memory or a USB key.

### Creating a Protocol


## B. Editing a Protocol

The **Edit** function allows changes to the denature temperature and time. Edited protocols may be saved under the existing name or be renamed (using the Save As function) to the controller memory or USB key.

1. Select Edit from the Home screen.
2. Choose the file location (controller or USB key) and the protocol to edit.
3. Make and confirm changes to the denature temperature and time
4. Save the edited file using the same name to memory or USB key.

To rename the file, select Save As and create a new name before saving to memory or a USB key.

### Editing a Protocol


## V. ONETEMP MODE

The **CytoBrite System** can be operated to maintain a single block temperature from 15 to 90°C using the OneTemp mode.

For information on creating and saving multi-temperature and multi-step protocols, see section **V. CREATING AND EDITING PROTOCOLS**.

1. From the Home screen, choose Run | OneTemp.
2. Set the desired block temperature from 15 to 90°C.
3. Press Select “Yes” to start the program. Set temperatures for the block will be displayed. Time will be listed as “Forever” as the program will continue until stopped by the user.
4. To end the program, use the **Forward** button to display the End Protocol screen and Select “Yes”. The instrument will display “OneTemp DONE” and total elapsed time.

### Running OneTemp Mode


## VI. RUNNING A PROTOCOL

The RUN function is used to run a previously created protocol stored in the memory or on a USB key.

1. Turn on the instrument. It preheats to 37C for five minutes before going to the main menu to run protocols etc.
2. Open the instrument lid and insert slide tray(s).
3. Select Run from the main menu.
4. Choose the file location (CytoBrite or USB key) and select the protocol.
5. The unit will rapidly heat the slides to the programmed denature temperature and time then rapidly cool to 37C.
6. Keep the slides in the unit at 37C for the desired incubation time. Alternatively, transfer the slide trays to a CytoBrite Slide Oven set at 37C.
7. End the incubate step by pressing any key to return to the main menu. The slide block remains at 37C.

### Running a Protocol

## VII. MANAGING FILES

### A. Viewing a Protocol


Saved protocols can be reviewed in read-only mode.

1. Select View from the Home screen.
2. Choose the file location (memory or USB) and protocol.
3. The programmed denature temperature and time is displayed.


### B. Copying a Protocol

A protocol can be copied and saved under a new name in the controller memory or USB key.

1. Select Files from the Home screen.
2. Choose the file location (controller or USB key) and protocol.
3. Select Copy.
4. Choose new file location and Select Yes.

### Copy/Rename/Delete Protocols

**C. Renaming a Protocol**

An existing file may be renamed and saved in the controller memory or USB key.

1. Select Files from the Home screen.
2. Choose the file location (controller memory or USB key).
3. Select Name.
4. Enter a new name by choosing characters (space; A-Z; 0-9) using the **Forward** button. Sequentially select a desired character using the **Select** button. Protocol names must be 8 characters in length, including spaces.
5. Select Yes to save.

**D. Deleting a Protocol**

1. Select Files from the Home screen.
2. Choose the file location (controller memory or USB key) and protocol.
3. Choose Delete and press **Select**.
4. Select Yes to confirm.

## VIII. SETTINGS AND PREFERENCES

### A. Datalogging to USB Key

A text file that captures the block time and temperature while a protocol is running can be created and written to the USB key.


1. From the Home screen, go to Settings | Datalog
2. Following the prompts on the screen, turn on datalogging and enter the data sampling frequency (5 sec to 60 min).
3. Insert a USB key into the USB DRIVE to capture a log file. While in use, the drive light will flash as data is transferred to the USB key.

A datalog file for a typical protocol is shown below:

```
Datalog for Protocol: 2MIN72C
Protocol Summary:
Denature 72C / 02mins 00secs
Incubate 37C / Forever
```

| Run Data Step | Time at Step | Set Temp | Actual Temp |
|---------------|--------------|----------|-------------|
| Denature      | 00:00:00     | 72C      | 36.7C       |
| Denature      | 00:00:15     | 72C      | 37.6C       |
| Denature      | 00:00:30     | 72C      | 43.1C       |
| Denature      | 00:00:45     | 72C      | 49.8C       |
| Denature      | 00:01:00     | 72C      | 56.3C       |
| Denature      | 00:01:15     | 72C      | 61.4C       |
| Denature      | 00:01:30     | 72C      | 66.1C       |
| Denature      | 00:01:45     | 72C      | 69.6C       |
| Denature      | 00:02:00     | 72C      | 71.3C       |
| Denature      | 00:02:15     | 72C      | 71.8C       |
| Denature      | 00:02:30     | 72C      | 71.9C       |
| Denature      | 00:02:45     | 72C      | 71.9C       |
| Denature      | 00:03:00     | 72C      | 71.9C       |
| Denature      | 00:03:15     | 72C      | 72.0C       |
| Denature      | 00:03:30     | 72C      | 72.1C       |
| Denature      | 00:03:45     | 72C      | 72.1C       |
| Incubate      | 00:04:00     | 37C      | 72.1C       |
| Incubate      | 00:04:15     | 37C      | 69.0C       |
| Incubate      | 00:04:30     | 37C      | 61.8C       |
| Incubate      | 00:04:45     | 37C      | 54.3C       |
| Incubate      | 00:05:00     | 37C      | 47.4C       |
| Incubate      | 00:05:15     | 37C      | 41.5C       |
| Incubate      | 00:05:30     | 37C      | 37.6C       |
| Incubate      | 00:05:45     | 37C      | 36.2C       |
| Incubate      | 00:06:00     | 37C      | 36.2C       |
| Incubate      | 00:06:15     | 37C      | 36.7C       |
| Incubate      | 00:06:38     | 37C      | 37.3C       |

### Creating a Datalog

### B. Calibrating Slide Block Temperature

The **CytoBrite System** comes calibrated from the factory to provide accurate slide temperatures from 15 to 90°C ± 0.5°C. Temperature accuracy should be checked periodically (according to your lab schedule) with a T-type NIST certified digital thermometer (SciGene cat. #1051-52-0).

The system is calibrated by entering offset values at 75°C and 37°C, as follows:

- Using the cable supplied with the T-type NIST certified digital thermometer (SciGene cat. # 1051-52-0), connect it to the blue jack on the rear of the instrument labeled “Temperature Calibration” and turn on power to the CytoBrite System.
- From the Home screen, go to the Settings menu and select “Calibrate”. The instrument will automatically heat to 75°C.
- When the temperature reaches 75°C on the controller display, wait 5 minutes.
- To determine the offset value, subtract the controller temperature from the thermometer temperature.
 


For example:

  - If the thermometer reads 75.5 and the controller reads 75.0, **add 0.5 to the existing offset value.**
  - If the thermometer reads 74.5 and the controller reads 75.0, **subtract 0.5 from the existing offset value.**
- Enter the new offset value in the Calibration screen by using the Forward or Back buttons on the controller, then press Select. The controller will now display temperature using new offset.
- Compare the new controller temperature to the thermometer.
 

If the temperatures match each other to within a few tenths of a degree, press Yes to save the new offset. If not, press No to enter a new offset value. Once the new offset is confirmed the instrument will automatically cool to 37°C.
- When the temperature reaches 37°C on the controller display, wait 5 minutes.
- Repeat steps 4 through 6 to determine a new offset value at 37°C.
- Confirm the second offset to return to the Settings menu.


**Calibration is complete.**

### Temperature Calibration

## C. Upgrading Software


The **CytoBrite System** displays the current software version on the controller when powering on the instrument. New software versions are installed via the USB DRIVE.

1. Turn on the instrument and insert the USB key with new software version into the USB Drive.
2. From the Main menu, go to Settings | Upgrade SW and press the Select button. A flashing progress indicator displays as the software is loading. Allow a minute for the upgrade to complete.
3. When finished, a message appears in the display to indicate success. Press Select to return to the Main Menu.

**The instrument is now ready to use.**

### Upgrading Software

## IX. MAINTENANCE

### A. Cleaning the Instrument

Clean outer metal surfaces and the heated lid using a soft cloth and mild, detergent-based cleaner. Avoid abrasive cleaners that can scratch surfaces.

Do **NOT** use caustic or strongly alkaline solutions (e.g., strong soaps, ammonia, or bleach at a higher concentration than specified above). If running radioactive or biohazardous reactions, consult your institution's radiation safety officer or biosafety office for assistance.

**Turn the power switch to the OFF position and unplug the power cord before performing any cleaning procedure!!**

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**Remove liquid or condensation (wipe dry) while instrument is OFF and block is cool or at room temperature!!**

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### B. Cleaning Air Vents

Clean the air intake (bottom of instrument) to prevent vents from becoming clogged. Remove any dust and debris, airflow may be hampered, effecting performance of the instrument and damage by overheating of components.

### C. Checking and Replacing Fuses


There are two fuses located in a removable fuse block below the power cord receptacle on the back of the unit. A blown fuse appears dark. Always replace fuses with those of the same amperage and voltage as shown on the label below the fuse block.

**Turn the power switch to the OFF position and unplug the power cord before performing any service procedure.**

To remove fuses:

1. Insert a small, flat blade screwdriver into the tab recess just below the plug receptacle.
2. Push down to release the fuse block.
3. Slide fuse holder out from the fuse block while holding the retaining tab out of the way.
4. Gently pry out the fuses.

## X. TROUBLESHOOTING

| Symptom                                | Cause                               | Solution   |
|--|-------------------------------------|--|
| Instrument does not turn on.           | Blown fuse(s)                       | Replace fuse(s) on back of unit, beneath power cord receptacle.                          |
| Block is not reaching set temperature. | Paper or other object blocking fans | Check underneath CytoBrite™ base for objects blocking the fans. Ensure proper clearance. |

## XI. SPECIFICATIONS

|                                      |  |
|--------------------------------------|--|
| <b>Electrical</b>                    |  |
| Cat. #2019-00-1                      | 115/230V AC; 50/60 Hz; 6A/3A           |
| <b>Dimensions</b> (cover closed)     |  |
| Outside (H x W x D)                  | 8 x 10 x 15 inches (20 x 25 x 38 cm)   |
| <b>Weight</b>                        |  |
| Instrument                           | 20 lbs (9 kg) net                      |
| <b>Performance</b>                   |  |
| Block Temperature Range              | 15 to 90°C                             |
| Block Temperature Regulation         | ± 0.2°C from set point                 |
| Block Temperature Accuracy           | ± 0.5°C                                |
| Slide-to-slide Temperature Variation | < 1.0°C                                |
| Heating Rate                         | 37 to 75°C in < 2 minutes              |
| Cooling Rate                         | 75 to 37°C in < 2 minutes              |
| <b>Controller</b>                    |  |
| Controller Type                      | Digital PID, single loop               |
| Display                              | Four line LED                          |
| Calibration Output                   | Thermocouple                           |
| <b>Environmental</b>                 |  |
| Ambient temperature                  | Operation: +15 to +32°C (58 to 90°F)   |
|                                      | Storage: -20 to +60°C (-4 to 140°F)    |
| Relative humidity                    | Operation: 20 to 80% RH non-condensing |
|                                      | Storage: 10 to 90% RH non-condensing   |
| Overvoltage category                 | II                                     |

## XII. ORDERING INFORMATION

| Cat. #    | Description  | UoM |
|-----------|--|-----|
| 2019-70-1 | CytoBrite Slide Oven, 115V. Includes ten 6-slide racks.                      | EA  |
| 2019-20-1 | CytoBrite Duo Slide Incubation System, 115/230V. Includes two 6-slide racks. | EA  |
| 2019-10-0 | CytoBrite Slide tray. Holds 1 to 6 slides.                                   | EA  |
| 2010-00-1 | FISH Wash Buffer 1 (0.4xSSC/0.3% IGEPAL, pH 7).                              | EA  |
| 2010-00-2 | FISH Wash Buffer 2 (2xSSC/0.1% IGEPAL, pH 7).                                | EA  |
| 2020-00-1 | CytoBond Removable Coverslip Sealant, 100 ml.                                | EA  |
| 2022-00-2 | CytoZyme Stabilized Pepsin, 50X Concentrate, 20 ml.                          | EA  |
| 2022-00-3 | CytoZyme Stabilized Pepsin, 50X Concentrate, 100 ml.                         | EA  |
| 2022-10-2 | CytoZyme Reaction Buffer, 1L.  | EA  |
| 2022-10-3 | CytoZyme Reaction Buffer, 4L.  | EA  |
| 2030-00-1 | Sodium Thiocyanate Pretreatment Reagent, 1L                                  | EA  |
| 2030-00-2 | Sodium Thiocyanate Pretreatment Reagent, 4L                                  | EA  |
| 1051-52-0 | Digital thermometer. Includes cable and NIST certificate                     | EA  |

## XIII. DECLARATION OF CONFORMITY

### CytoBrite Slide Incubation System

SciGene  
 1287 Reamwood Ave  
 Sunnyvale, CA 94089 USA

Declares that the above referenced product(s) meets the essential requirements of the following European Union Directives by using the relevant standards shown below to indicate compliance.

|                                  |             |         |   |
|----------------------------------|-------------|---------|---|
| <b>EMC Directive 2004/108/EC</b> |             |         |   |
| <b>EN 61326-1</b>                | <b>2013</b> |         | <b>Electrical equipment for measurement, control and laboratory use to include:</b> |
| EN 55011                         | 2010        | Class A |   |
| EN 61000-3-2                     | 2009        |         |   |
| EN 61000-3-3                     | 2013        |         |   |
| EN 61000-4-2                     | 2009        |         |   |
| EN 61000-4-3                     | 2006        |         |   |
| EN 61000-4-4                     | 2012        |         |   |
| EN 61000-4-5                     | 2006        |         |   |
| EN 61000-4-6                     | 2014        |         |   |
| EN 61000-4-8                     | 2010        |         |   |
| EN 61000-4-11                    | 2004        |         |   |

|                                 |             |  |  |
|---------------------------------|-------------|--|--|
| <b>LVD Directive 2006/95/EC</b> |             |  |  |
| <b>IEC 61010-1</b>              | <b>2013</b> |  | <b>Safety requirements for measurement, control and laboratory use</b> |
|                                 |             |  | <b>3rd Edition Part 1: General requirements</b>                        |

### RoHS Directive 2011/65/EU

Terry Gill  
 Name of Authorized Representative

Sunnyvale, California, USA  
 Place of Issue

Director of Product Manufacturing  
 Title of Authorized Representative

March 19, 2015  
 Date of Issue

  
 Signature of Authorized Representative