# CytoBrite® Duo

## Slide Incubation System

### USER MANUAL

Cat. #2019-20-1 (115/230V)

FOR RESEARCH USE ONLY

---

**Serial Number**

The following serial number identifies the specific instrument you have purchased and must be referenced when requesting service. A copy is affixed to the instrument.

Technical Service: (408) 733-7337, [techserv@scigene.com](mailto:techserv@scigene.com)

**Warranty**

SciGene warrants that the heating unit described in this manual shall be free of defects in materials and workmanship for a period of 12 months from date of delivery. This warranty does not cover removable slide trays or accessories. In the event of a defect during the warranty period, SciGene's limit of liability will be to provide replacement parts at no charge or, at its sole discretion, replace the product. The foregoing warranty is void in the event the unit was abused or modified or used in a manner inconsistent with its intended purpose. SciGene makes no other warranty, expressed or implied including warranties of merchantability and fitness for a particular purpose. In no event shall SciGene be liable for any direct, indirect, special, incidental or consequential damages or for any damages resulting from loss arising out of or in connection with the sale, use or performance of the product.

**Copyright**

Copyright ©2014-2015 SciGene Corporation. All rights reserved. SciGene, CytoBrite and CytoBond are trademarks of SciGene Corporation, Sunnyvale, CA. All other trademarks used in this manual are the property of their respective owners.

## Table of Contents

<b>I. SAFETY NOTICES</b> .....	<b>3</b>
A. Intended Use .....	3
B. Instrument Safety .....	3
C. Symbols and Conventions .....	3
D. Warnings .....	4
E. Cautions .....	5
F. Electromagnetic Compatibility (EMC) .....	5
G. FCC Warnings and Notes .....	5
H. Compliance .....	6
<b>II. UNPACKING AND SET UP</b> .....	<b>7</b>
A. Unpacking the System .....	7
B. Items Provided .....	7
C. Coverslip Sealant Requirements .....	7
D. Environmental Requirements .....	7
E. Voltage Selection and Fuse Type .....	7
F. Lid Configuration .....	8
G. Installation .....	8
<b>III. USING YOUR CYTOBRITE DUO SYSTEM</b> .....	<b>9</b>
A. Main Components, Controls and Accessories .....	9
B. Voltage Selector and Connections .....	10
C. Using CytoBond Removable Coverslip Sealant .....	10
D. Handling Slide Trays .....	10
E. Using the Temperature Controllers .....	11
F. Calibrating the Temperature Controllers .....	11
<b>IV. MAINTAINING YOUR CYTOBRITE DUO SYSTEM</b> .....	<b>12</b>
A. Powering Off .....	12
B. Cleaning the Instrument .....	12
C. Checking and Replacing Fuses .....	12
<b>V. TROUBLESHOOTING</b> .....	<b>13</b>
<b>VI. SPECIFICATIONS</b> .....	<b>13</b>
<b>VII. RELATED INSTRUMENTS, REAGENTS AND ACCESSORIES</b> .....	<b>13</b>
<b>VIII. DECLARATION OF CONFORMITY</b> .....	<b>14</b>

## I. SAFETY NOTICES

### A. Intended Use

The **CytoBrite® Duo Slide Incubation System** is a dual temperature slide heating instrument for FISH/ISH applications. The instrument should only be used according to the instructions provided in this manual. If the equipment is used in a manner not specified by the manufacturer (SciGene), protection provided by the equipment may be impaired.

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### B. Instrument Safety

Before operating the instrument, read the information in this section concerning hazards and potential hazards. Ensure that anyone involved with the instrument's operation is instructed in both general safety practices for laboratories and specific safety practices for the instrument.

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### C. Symbols and Conventions

The following chart is an illustrated glossary of the electrical symbols that are used on the **CytoBrite Duo Slide Incubation System**. Whenever such symbols appear on instruments, please observe appropriate safety measures.

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#### 1. Electrical Symbols

This symbol indicates that this is a protected ground terminal that must be connected to earth ground before any other electrical connections are made to the instrument.

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CAUTION: This symbol alerts you to consult this User Manual for further information and to proceed with caution.

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This symbol indicates the OFF position of the main POWER switch.

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This symbol indicates the ON position of the main POWER switch.

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## 2. Non-Electrical Symbols

CAUTION: This symbol illustrates a heat hazard. Proceed with caution when working around these areas to avoid being burned by hot components.

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CAUTION: This symbol alerts you to consult this Operator's Manual for further information and to proceed with caution.

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## D. Warnings

Failure to comply with the following warnings that are affixed to the product can lead to possible personal injury or death.

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This symbol on the rear of the instrument indicates the presence of the fuse box. **Warning: For Continued Protection Against Fire, Replace Only with Same Type Rating of Fuse.** Always disconnect the power cord before attempting to replace the fuse.

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## E. Cautions

Failure to comply with the following cautionary statement affixed to the product may lead to possible personal injury.

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This symbol located in front of the heat blocks indicates the potential presence of a Hot Surface. Use care when working in this area to avoid being burned.

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This symbol located on the back of the instrument warns the user to select the correct voltage before connecting the power cord. Operating with the voltage selector switch in the wrong position will damage the instrument and create a safety hazard.

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## F. Electromagnetic Compatibility (EMC)

The **CytoBrite Duo Slide Incubation System** is a Class A digital device under FCC Title 47 Part 15B and designated as Class A electrical equipment for measurement, control, and laboratory use (EN61326).

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## G. FCC Warnings and Notes

This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. The instrument qualifies as an “exempted device” under 47 CFR 15.103(c), in regard to the cited FCC regulations in effect at the time of manufacture. Shielded cables must be used with this unit to ensure compliance with the Class A FCC limits.

## H. Compliance

### 1. Electromagnetic Compatibility (CE)

All instruments shipped to the European Union (EU; formerly known as the European Community) have the “CE” label on the back of the instrument, signifying that these instruments comply with the Electromagnetic Compatibility and Low Voltage Directives.

### 2. Restriction of Hazardous Substances (RoHS)

The RoHS directive was adopted in February 2003 by the European Union. It restricts (with exceptions) the use of six hazardous materials (including lead and mercury) in the manufacture of various types of electronic equipment. All SciGene instruments carrying the CE mark are also RoHS compliant.

### 3. Electrical and Electronic Equipment Waste (WEEE)

WEEE is a European Community directive that became European Law in February 2003. The crossed out wheeled bin symbol on a product or its packaging indicates that it must NOT be disposed of with standard waste at the end of its lifetime. Instead, it is your responsibility to return it to a designated collection point for electronic equipment recycling. Contact your local SciGene distributor for information on drop off locations.

**II. UNPACKING AND SET UP**

**A. Unpacking the System**

The **CytoBrite Duo Slide Incubation System** is shipped in a single carton with empty slide trays placed in the heat blocks. Remove any packing material around the instrument and lift out; taking care not to damage the hinged lids. Carefully inspect for damage.

*If the instrument is damaged upon arrival, please retain all shipping materials and contact SciGene or your local distributor for assistance.*

**B. Items Provided**

The following items are included with the CytoBrite Duo instrument:

- 2x Slide Trays
- Power Cord
- User Manual

**C. Coverslip Sealant Requirements**

Ready-to-use **CytoBond Removable Coverslip Sealant** (SciGene cat. #2020-00-1) is required for use of the **CytoBrite Duo System**. **CytoBond Sealant** is used to temporarily seal coverslips to slides, providing an evaporative seal over prolonged incubations and high temperatures without humidification.

**D. Environmental Requirements**

Ensure that the area where the system is installed meets the following conditions, for reasons of safety and performance:

Ambient temperature	+15 to +32°C (58 to 90°F)
Relative humidity	20 to 80% RH non-condensing
Air flow clearance	3 inches (8 cm) minimum on both sides

**E. Voltage Selection and Fuse Type**


The **CytoBrite Duo System** may be used with either 115V (US, Canada, Japan) or 230V electricity by sliding the voltage selector on the back of the instrument. Operating the system with the switch in the wrong position may damage the instrument and create a safety hazard. Before connecting the power cord verify that the switch is set to the correct voltage for your locality with the correct fuses installed:

- 115 Voltage requires two 4 Amp fuses
- 230 Voltage requires two 2 Amp fuses

See section IV. C. **Checking and Replacing Fuses** for fuse replacement instructions.

**To avoid a safety hazard or damage to the instrument, select the correct voltage and ensure the correct fuses are installed for your locality BEFORE connecting the power cord.**


**F. Lid Configuration**

The lids covering the slide trays are configured at the factory to open to the right but can be easily switched to open in the opposite direction.


Follow these steps to re-configure one or both lids to open to the left:

1. Remove the bottom two screws holding the lid to the bracket on the right side of the instrument. Carefully lift the lid and set aside.
2. Remove the corresponding two screws on the left side of the instrument to release a small cover plate.
3. Turning the lid around, place it back on the instrument and secure it by replacing the bottom two screws of the bracket.
4. On the right side of the instrument, re-install the small cover plate with the two screws previously removed. Repeat steps 1 through 4 for the second lid, if desired.

**G. Installation**


Place the instrument on a level surface within a few feet of the power source. Position the system so the lids can be easily opened without interference. Allow 18 inches of vertical clearance and 16 inches of horizontal clearance. Ensure 3 inches of clearance at the back to attach the power cord. Using only the power cord provided, plug it into the back of the unit and then to a properly grounded outlet. Turn on power to the block heaters using the ON/OFF switches on the front control panel.

III. USING YOUR CYTOBRITE DUO SYSTEM


A. Main Components, Controls and Accessories

Name	Function
<b>Power Switches</b>	Turns ON/OFF power to the instrument.
<b>Temperature Controllers</b>	Used to set, control and observe block temperatures.
<b>Hinged Lids</b>	Covers the slides during sample incubation.
<b>Heat Blocks (front and back)</b>	Heats slides held in trays.
<b>Slide Trays (2 included)</b>	Holds slides for incubation.
<b>CytoBond® Removable Coverslip Sealant</b>	Seals coverslips to slides without humidification.
<b>Digital Thermometer (optional)</b>	Used to verify and calibrate block temperatures.

**B. Voltage Selector and Connections**

Name	Function
<b>Calibration Jacks</b>	Connects T-type thermometer to internal temperature sensor
<b>Voltage Selector</b>	Used to switch between 115 and 230 Volts for your locality
<b>Power Connection</b>	Connection for provided power cord
<b>Fuse Block Holder</b>	Holds removable block for replacing fuses

**C. Using CytoBond Removable Coverslip Sealant**


The **CytoBrite Duo System** eliminates the need for humidification with the use of **CytoBond Removable Coverslip Sealant** (SciGene cat. #2020-00-1). CytoBond is a gel-like adhesive applied to coverslips that temporarily seals them to slides.

Follow the steps below to prepare slides with **CytoBond Sealant**:

1. Apply a continuous bead around coverslips.
2. Wait 5 to 30 minutes.
3. Denature/incubate using your times/temperatures.
4. Remove in one clean piece with forceps.

**D. Handling Slide Trays**


Two slide trays are included with the **CytoBrite Duo System**. Each removable tray holds 1 to 6 slides from assay setup to hybridization. To load a tray, hold its end tabs only (to avoid touching slides) and place the tray firmly onto the desired heat block such that all slides remain flat and level, in full contact of the block. Failure to properly insert the slide tray or slides may adversely affect performance.

Each tray holds 1 to 6 slides

### E. Using the Temperature Controllers

The CytoBrite Duo System is equipped with two Watlow EZ-Zone temperature controllers. Each controller has two LED displays and four push buttons. The left display shows ACTUAL block temperature in red. The right display shows the SET or programmed temperature in green. To enter the SET temperature, simply push the up and down arrows until the desired temperature is shown in the green display on the right. The unit will adjust the heat of the block until the SET temperature is attained.

### F. Calibrating the Temperature Controllers

The controllers are calibrated at the factory. The front block is calibrated at 37°C and the back block at 75°C. Calibration is required ONLY if the block temperature differs by more than one degree (1°C) from the actual temperature on the controller when using a NIST-calibrated T-type digital thermometer (sold separately, SciGene cat. #1051-52-0).


Follow these steps to adjust a controller to achieve accurate temperatures:

1. For each block, turn ON the controller and set to the desired temperature. Allow 30 minutes to stabilize.
2. Using the thermometer cable provided, plug one end into the corresponding blue jack on the back of the system and the other into the digital thermometer.
3. Turn on the thermometer and allow 1 minute to stabilize. Actual temperature is displayed.
4. Calculate the difference between the thermometer and the controller to determine the adjustment value. For example, if the controller is set to 75.0°C and the thermometer reads 73.9°C while the controller reads 75.0°C, then the adjustment value is -1.1°C.
5. On the controller, press the up and down arrows simultaneously for 3 seconds. The left display shows "A1" and the right display shows "open".
6. Press the Advance Button (green circle) 3 times until the right display shows "i.CA". The left display will show the offset value when the unit was last calibrated.
7. Using the up or down arrows, add the calculated value from step 4 to the existing offset value. For example, if the adjustment value is -1.1°C and the current offset is -0.3°C then the new offset is -1.4°C.
8. Press the Infinity Key ( $\infty$ ) twice to exit calibration and return to the operation display. Verify that the thermometer matches the controller. Your system is now calibrated to provide accurate temperatures between ambient +5 and 90°C.

## IV. MAINTAINING YOUR CYTOBRITE DUO SYSTEM

### A. Powering Off

Turn OFF power to both front and back block heaters and unplug the power cord before performing any service procedure.

### B. Cleaning the Instrument

Clean outer metal surfaces and the heated lid using a soft cloth and mild, detergent-based cleaner. Avoid abrasive cleaners that can scratch surfaces. Do **NOT** use caustic or strongly alkaline solutions (e.g., strong soaps, ammonia, or bleach). Clean the heat blocks to avoid buildup of residue that can interfere with the heat control system.

If running radioactive or biohazardous reactions, consult your institution's radiation safety officer or biosafety office for assistance.

**Turn both power switches to the OFF position and unplug the power cord before performing any cleaning procedure!**

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**Remove liquid or condensation (wipe dry) while instrument is OFF and at room temperature!**

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### C. Checking and Replacing Fuses

Two fuses are located in the removable block below the power cord receptacle on the back of the unit. Blown fuses will have a broken filament or appear dark. Always replace fuses with the same amperage and voltage as shown on the label below the fuse block.


**Turn both power switches to the OFF position and unplug the power cord before performing any service procedure.**

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To remove fuses:

1. Insert a small, flat blade screwdriver into the tab recess just below the plug receptacle.
2. Push down to release the fuse block.
3. Slide fuse holder out from the fuse block while holding the retaining tab out of the way.
4. Gently pry out the fuses.

## V. TROUBLESHOOTING

Symptom	Likely Cause	Solution
Instrument does not turn on.	Blown fuse(s)	Replace fuse(s) on back of unit, beneath power cord receptacle.
Controller is not responding when buttons are pressed.	Controller is faulty or is not properly programmed	Cycle unit on and off. Contact SciGene or your local distributor for support.

## VI. SPECIFICATIONS

Electrical	
Cat. #2019-20-1	115/230V AC; 50/60 Hz; 4A/2A
Dimensions and Weight	
Outside (lids closed) H x W x D	6 x 10 x 13.5 inches (15 x 25.5 x 34 cm)
Weight	12 lbs (5.4 kg)
Performance and Controls	
Temperature Range	Ambient +5°C to 90°C
Temperature Regulation	± 0.2°C
Slide-to-slide Temp Variation	< 0.5°C
Preheat Time	< 15 minutes
Temperature Controllers	Digital PID
Temperature Displays	Actual or Set single LEDs
Digital Thermometer Output	T-Type Thermocouple

## VII. RELATED INSTRUMENTS, REAGENTS AND ACCESSORIES

Cat. #	Description	UoM
2019-00-1	CytoBrite Slide Incubation System, 115V/230V. Includes two 6-slide racks.	EA
2019-70-1	CytoBrite Slide Oven, 115V. Includes ten 6-slide racks.	EA
2019-10-0	CytoBrite Slide tray. Holds 1 to 6 slides.	EA
2010-00-1	FISH Wash Buffer 1 (0.4xSSC/0.3% IGEPAL, pH 7).	EA
2010-00-2	FISH Wash Buffer 2 (2xSSC/0.1% IGEPAL, pH 7).	EA
2020-00-1	CytoBond Removable Coverslip Sealant, 100 ml.	EA
2022-00-2	CytoZyme Stabilized Pepsin, 50X Concentrate, 20 ml.	EA
2022-00-3	CytoZyme Stabilized Pepsin, 50X Concentrate, 100 ml.	EA
2022-10-2	CytoZyme Reaction Buffer, 1L.	EA
2022-10-3	CytoZyme Reaction Buffer, 4L.	EA
2030-00-1	Sodium Thiocyanate Pretreatment Reagent, 1L	EA
2030-00-2	Sodium Thiocyanate Pretreatment Reagent, 4L	EA
1051-52-0	Digital thermometer, T-type. Includes cable and NIST certificate	EA

**VIII. DECLARATION OF CONFORMITY**

**CytoBrite® Duo Slide Incubation System**

SciGene  
 1287 Reamwood Avenue  
 Sunnyvale, CA 94089 USA

Declares that the above referenced product(s) meets the essential requirements of the following European Union Directives by using the relevant standards shown below to indicate compliance.

**EMC Directive 2004/108/EC**

**EN 61326-1      2013      Electrical equipment for measurement, control and laboratory use to include:**

EN 55011	2010	Class A
EN 61000-3-2	2009	
EN 61000-3-3	2013	
EN 61000-4-2	2009	
EN 61000-4-3	2006	
EN 61000-4-4	2012	
EN 61000-4-5	2006	
EN 61000-4-6	2014	
EN 61000-4-8	2010	
EN 61000-4-11	2004	

**LVD Directive 2006/95/EC**

**EN 61010-1      2013      Safety requirements for measurement, control and laboratory use  
 3<sup>rd</sup> Edition Part 1: General requirements**

**RoHS Directive 2011/65/EU**

Terry Gill  
 \_\_\_\_\_  
 Name of Authorized Representative

Sunnyvale, California, USA  
 \_\_\_\_\_  
 Place of Issue

Director of Product  
 Manufacturing  
 \_\_\_\_\_  
 Title of Authorized Representative

March 13, 2015  
 \_\_\_\_\_  
 Date of Issue

  
 \_\_\_\_\_  
 Signature of Authorized Representative